





Reforma

ExtractionBasket



Pancreatic Stone Extraction Basket “Reforma”


Reforma is used for extraction of pancreatic stone and bile duct stone with Endoscopic procedure.



Wire O.D.:0.05mm
consist of 72 wires

- **Basket part made of shape memory alloy is suitable for internal wall of duct. Reforma catch the stone, (& sludge) in side pancreatic duct and bile duct.**
- **Control the handle to deploy basket and possible to manipulate opening and closing basket.**

Product Feature



Catalog No.	Basket Diameter	Number of Basket Wire	Outer Diameter of Sheath Tip	Outer Diameter of Sheath hand	Guidewire Used with	Length
BF72-05E	5mm	72	1.90mm (5.7Fr)	2.45mm (7.4Fr)	0.025inch	1910mm