

Trefle

Newly designed biliary cytology brush for biopsy



PIOLAX
Endoscopy

Biliary Cytology brush for Biopsy “Trefle”


Approach of New Device for diagnosis for Bile Duct Stricture

Intended for use


Trefle harvests the cells and tissue fragment off from Stenosis site of Pancreatic duct (bile duct) for using diagnosis.

3 Key Points of “Trefle”

- Use by Guide wire guided
- Scratching by not Brush, metallic loops.
- Not only brushing tissues and collecting tissues **can** be performed for tissue diagnosis.



Product Feature



Catalog No.	Sheath outer Diameter	Length	Maximum Diameter of Guidewire	Number of loop in Biopsy	Deployment Diameter of Loop in Biopsy
PB7-3L5S	2.45mm (7.3Fr)	1950mm	0.89mm (0.035inch)	3	4.5mm