

# Uneven Double Lumen Cannula

The Expansion for new ERCP procedure by Double Lumen Cannula





## ERCP Cannula “Uneven Double Lumen Cannula”


### Unique Distal Tip



### Multi Usage



### Long nose



### Short nose

During using G.W., can perform injection and suction(Cytotechnology) .

2<sup>nd</sup> G.W. replacement or exchanging G.W. (0.025⇔0.035) can be performed.

| O.D.<br>(Distal Tip) | Nose<br>Length | O.D.<br>(Proximal Tip) | Length | Compatible Guide Wire         |
|----------------------|----------------|------------------------|--------|-------------------------------|
| 3.6Fr                | 5mm            | 6Fr.                   | 1830mm | Distal Lumen:0.66mm(0.025")   |
|                      | 30mm           |                        |        | Proximal Lumen:0.80mm(0.035") |

